




DONJOY ULTRASLING® III PATIENT APPLICATION

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO PROPER DEVICE FUNCTIONING.



INTENDED USE AND INDICATIONS

Postoperative rotator cuff repairs, Bankart procedures, capsular shifts, glenohumeral dislocation/subluxation, fractures, soft tissue strains/repairs, and joint replacements.

APPLICATION INFORMATION

1. Detach shoulder strap buckles and open front panel. Position elbow in the sling as far back as possible.
 - a. To maintain proper positioning, elbow should be at 90 degrees of flexion and pillow should rest at side of body (DO NOT have pillow rest on stomach).
2. Place the shoulder strap over the opposite shoulder. Connect the shoulder strap to the sling using the two quick release buckles. For exercise and rehab detach the outside quick release buckle to increase the range of motion.
3. Secure strap at the top of the sling. If the thumb strap is desired, attach strap at the front of the sling (thumb should rest on strap comfortably without added pressure).
4. Ensure the pillow is at the waistline of the affected side and buckle the waist strap to the pillow.
5. Follow physician's orders for duration and frequency of use.
6. To remove: unbuckle waistband quick release buckle at front of pillow. Un-Velcro top closure strap and thumb strap from sling. Unbuckle two quick release buckles on shoulder strap. *Do not Un-Velcro shoulder strap or waistband unless sizing needs to be adjusted. Do not remove sling from pillow unless instructed by physician.*

WARNINGS & PRECAUTIONS

If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, consult your medical professional immediately. Never wear the support over an open wound.

CLEANING INSTRUCTION

Wipe clean with a wet washcloth and mild detergent or baby shampoo. Wipe again with wet cloth and AIR DRY only; do not heat dry. Regular cleaning is recommended.

Midwest Orthopaedics at Rush DME Store Contact Information

Please call or stop in if you have any problems with sling fit or if you have any other questions

Chicago - 1611 W Harrison St. Suite 103, Chicago, IL - (312) 432-2482

Westchester - 2450 S Wolf Road Suite G, Westchester, IL - (708) 273-8426

Oak Park - 610 S Maple Ave. Suite 1550, Oak Park, IL - (312) 432-2550

Munster - 9200 Calumet Ave. Suite S300, Munster, IN - (708) 492-5450

Oak Brook - 2011 York Road, Oak Brook, IL - (708) 492-5664

Naperville - 55 Shuman Blvd Suite 178, Naperville, IL - (630) 339-2262

Joliet - 963 129th Infantry Dr. Suite 100, Joliet, IL - (708) 492-5921